



Recombinant Chicken Hyaluronan synthase 3 (HAS3)

Product Code	CSB-CF010142CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O57425
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	SYFGCVQCISGPLGMYRNTLLQQFLEDWYHQTFGLGKCSFGDDRHLTNRVLS LG YQTKYT PRSRCLTETPTRYLRWLNQQTRWSKSYFREWLYNALWFHKHHLWMTYESVV TGFFPFFLI ATVIQLFYRGRINILLFLLTVQLVGIIKATYACFLRGS AEMIFVSLYALLYMSSLL PAK MFAIATINKS
Source	in vitro E.coli expression system
Target Names	HAS3
Protein Names	Recommended name: Hyaluronan synthase 3 EC= 2.4.1.212 Alternative name(s): Hyaluronate synthase 3 Hyaluronic acid synthase 3 Short name= CHAS3 Short name= HA synthase 3
Expression Region	1-190
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.