



Recombinant Zea mays Cell number regulator 1 (CNR1)

Product Code	CSB-CF469624ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B6TZ45
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	MYPSAPPDAYNKFSAGAPPTAPPPPAAYHQQQQHGANMDTSRPGGGLRK WSTGLFHCM DPGNCLITCLCPCVTFGQVADIVDKGTCPCIASGLVYGLICASTGMGCLYSCLY RSKLRA EYDLDEGECPDILVHCCCEHLALCQEYRELKNRGFDLGIGWEANMDRQRRGV AGGGAVMG APPAIPLGMIR
Source	in vitro E.coli expression system
Target Names	CNR1
Protein Names	Recommended name: Cell number regulator 1 Alternative name(s): ZmCNR01
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.