



Recombinant Human herpesvirus 6B Protein U15 (U15)

Product Code	CSB-CF868743HKA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9QJ48
Product Type	Transmembrane Protein
Immunogen Species	Human herpesvirus 6B (strain Z29) (HHV-6 variant B) (Human B lymphotropic virus)
Sequence	MDVWKRQRLQECRELCPLPALMSLSNILSNTTEIIVKYLKMFSTMYRFILPAL TLSMT VTKSVVIEMLFILKRWEEINQFFRLNIRKVND CFVVAQFTNIPVKRKIIVLLYMLTS RQE KQLFLNMIYAFLEKSHLRLGDDEEQNAIRFFSYVDELHLTRDVLLEIIVKLNTEI NQTM ELLLSYNELAG
Source	in vitro E.coli expression system
Target Names	U15
Protein Names	Recommended name: Protein U15
Expression Region	1-191
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.