



Recombinant TM2 domain-containing protein C02F5.13 (C02F5.13)

Product Code	CSB-CF351020CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P61228
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	NNEFRIEFYPNNEKSEKCFDSSKENDLLDLFYVSTNPLGPVVECRFLENSFIL CEDPVP LYGPGQTGQQPANESFRNEGKCLKMGGYRAEDVEFTNVKCRVLPICIECHGP RTFTKSTPC IYNGHYFLTLLYSIFLGVVAVDRFCLGYSAMAVGKLMTLGGFGIWWIVDIFLLV LGVL GPADDSSWEPYY
Source	in vitro E.coli expression system
Target Names	C02F5.13
Protein Names	Recommended name: TM2 domain-containing protein C02F5.13
Expression Region	19-210
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.