



Recombinant Rat Apoptosis regulator BAX (Bax)

Product Code	CSB-CF720536RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q63690
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MDGSGEQLGGGGPTSSEQIMKTGAFLLQGFIQDRAGRMAGETPELTLEQPPQ DASTKKLS ECLRRIGDELDSNMELQRMIADVDTDSPREFFRVAADMFA DGNFNWGRVVA LFYFASKL VLKALCTKVPELIRTIMGWTLDFLRERLLVWIQDQGGWDGLLSYFGTPTWQTV TIFVAGV LTASLTIWKKMG
Source	in vitro E.coli expression system
Target Names	Bax
Protein Names	Recommended name: Apoptosis regulator BAX
Expression Region	1-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.