



Recombinant Human Transmembrane protein 11, mitochondrial (TMEM11)

Product Code	CSB-CF023679HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P17152
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAAWGRRRLGPGSSGGSARERVLSATDCYIVHEIYNGENAQDQFEYELEQA LEAQYKYI VIEPTRIGDETARWITVGNCLHKTAVLAGTACLFPLALPLDYSHYISLPAGVLS LACCT LYGISWQFDPCKYQVEYDAYKLSRLPLHTLTSSTPVVLRKDDLHRKRLHNTI ALAALV YCVKKIYELYAV
Source	in vitro E.coli expression system
Target Names	TMEM11
Protein Names	Recommended name: Transmembrane protein 11, mitochondrial Alternative name(s): Protein PM1 Protein PMI
Expression Region	1-192
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.