



Recombinant *Drosophila yakuba* Calcium channel flower (flower)

Product Code	CSB-CF455787DMR
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	B4PD01
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila yakuba</i> (Fruit fly)
Sequence	MSFAEKITGLLARPNQQDPIGPEQPWYLKYGSRLLGIVAAFFAILFGLWNVFSII TLSVS CLVAGIIQMVAGFVVMLLEAPCCFVCFEQVNVIAADKVDKPLYFRAGLYITMAIP PIILC FGLASLFGSGLIFGTGVVYGMMALGKKASAEDMRAAAQQTFGGNTPAQTNDR AGIVNNAQ PFSFTGAVGTDSNV
Source	in vitro E.coli expression system
Target Names	flower
Protein Names	Recommended name: Calcium channel flower
Expression Region	1-194
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.