



Recombinant *Drosophila melanogaster* Mitochondrial import inner membrane translocase subunit Tim22 (Tim22)

Product Code	CSB-CF812853DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8IN78
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MSVLPNAGNNESSVTAPKMFGDPDLDRMAMEYVGNMQRYRENIIPKSNGPV KIKTNEEK FIETAVESCGFKCTMACVMGYGLGAALGLFSASVNPNMADPFANEKKQTARE VFREMRST THSYAKNFALIGCVFSAVECTIESHRGVTDWKNNGTYAGGITGGLIGLRAGVKA GIIGGLG FAAFSTADYYMYSR
Source	in vitro <i>E.coli</i> expression system
Target Names	Tim22
Protein Names	Recommended name: Mitochondrial import inner membrane translocase subunit Tim22
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.