



Recombinant Human Neurensin-1 (NRSN1)

Product Code	CSB-CF815587HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8IZ57
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MSSCSNVCGSRQAQAAAEGGYQRYGVRSYLHQFYEDCTASIWEYEDDFQIQ RSPNRWSSV FWKVGLISGTVFVILGLTVLAVGFLVPPKIEAFGEADFVVVDTHAVQFNSALDM YKLAGA VLFCIGGTSMAGCLLMSVFKSYSKEEKFLQQKFKERIADIKHAHTQPVTKAPGP GETKIP VTLSRVQNVQPLLAT
Source	in vitro E.coli expression system
Target Names	NRSN1
Protein Names	Recommended name: Neurensin-1 Alternative name(s): Neuro-p24 Vesicular membrane protein of 24 kDa Short name= Vesicular membrane protein p24
Expression Region	1-195
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.