



Recombinant *Hylobates syndactylus* Cytochrome b-c1 complex subunit Rieske, mitochondrial (UQCRFS1)

Product Code	CSB-CF731660HXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q69BJ6
Product Type	Transmembrane Protein
Immunogen Species	<i>Symphalangus syndactylus</i> (Siamang) (<i>Hylobates syndactylus</i>)
Sequence	SHTDVKVPDFSEYRRPEVLDSTKSSRESSEARKGFSYLVTAVTTVGVAYAAKN AVTQFVS SMSASADVLALAKIEIKLSDIPEGKNMAFKWRGKPLFVRHRTQKEIEQEAAVEL SQLRDP QHDLDRVKRPEWVILIGVCTHLGCVPIANAGDFGGYYCPCHGSHYDASGRIRL GPAPLNL EVPTYEFTSDDMVIVG
Source	in vitro E.coli expression system
Target Names	UQCRFS1
Protein Names	Recommended name: Cytochrome b-c1 complex subunit Rieske, mitochondrial EC= 1.10.2.2 Alternative name(s): Complex III subunit 5 Cytochrome b-c1 complex subunit 5 Rieske iron-sulfur protein Short name= RISP Ubiquinol-cytoc
Expression Region	79-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.