



Recombinant Gorilla gorilla gorilla Cytochrome b-c1 complex subunit Rieske, mitochondrial (UQCRFS1)

Product Code	CSB-CF731661GGZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q69BK4
Product Type	Transmembrane Protein
Immunogen Species	Gorilla gorilla gorilla (Western lowland gorilla)
Sequence	SHTDIKVPDFSEYRRLEVLDSTKSSRESSEARKGFSYLVTVGVTTVGVAYAAKN AVTQFVS SMSASADVLALAKIEIKLSDIPEGKNMAFKWRGKPLFVRHRTQKEIEQEAAVEL SQLRDP QHDLDRVKKPEWVILIGVCTHLGCVPIANAGDFGGYYCPCHGSHYDASGRIRL GPAPLNL EVPTYEFTSDDMVIVG
Source	in vitro E.coli expression system
Target Names	UQCRFS1
Protein Names	Recommended name: Cytochrome b-c1 complex subunit Rieske, mitochondrial EC= 1.10.2.2 Alternative name(s): Complex III subunit 5 Cytochrome b-c1 complex subunit 5 Rieske iron-sulfur protein Short name= RISP Ubiquinol-cytoc
Expression Region	79-274
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.