



Recombinant Danio rerio Mpv17-like protein 2 (mpv17l2)

Product Code	CSB-CF728529DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6DGV7
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MISRQGKEFLVRIAGYWKPLFRGRFLIVTNTVSCGGMLAAGDLIQQTREIRRTP GRTRDW SRTGCMFAVGCSMGPFMHWYQWLDKYFIGNGINNVCKKVLVDQLVASPTL GAWYFLGMG MMEGHTFIEAQQEFRDKFWEFYKADWCVWPAAQMINFYFLPPKFRVLYVNIV TLGWDTYL SYLKHRDTVEVTKEADGTA
Source	in vitro E.coli expression system
Target Names	mpv17l2
Protein Names	Recommended name: Mpv17-like protein 2
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.