



Recombinant Drosophila melanogaster Signal peptidase complex subunit 2 (Spase25)

Product Code	CSB-CF895680DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VYY2
Product Type	Transmembrane Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Sequence	MGKKDEKSQQGEELVKVNWGSAVKHALDDAVKTCLLGDRPQLKEQFGLV NTRLALCAL AVSVAIMAHAWDFTHPFPE SRPVLLFSVLAYFALLGILTLHSSFREKGTFAVAL QDKER ERLWEASSDMRKYDDKYLLTSLVRDTKNGKRREQSSNKSCAAFIDQNGIVLD NLVANEVN RLFNALAADKKNASSLSSN
Source	in vitro E.coli expression system
Target Names	Spase25
Protein Names	Recommended name: Signal peptidase complex subunit 2 EC= 3.4.-.- Alternative name(s): Microsomal signal peptidase 25 kDa subunit Short name= SPase 25 kDa subunit
Expression Region	1-199
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.