



Recombinant Schizosaccharomyces pombe Mitochondrial import inner membrane translocase subunit tim21 (tim21)

Product Code	CSB-CF527241SXV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O94618
Product Type	Transmembrane Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Sequence	SLASDSLNKQRKPQEEGLRLARIFKDPSPNKAWKDLTAPQKAYRTSANIGNFSIVI FGGGVF GLIYALVTSIWKGEAHYGDEAFELLKANEECRYVFGDHMKALGEATHPLRRTH GILTSR VWDHHGVEHLVLQFHLIGNERKGHVFGRLVNVQGDYKWEYLFVDVANYGKIII FDHTNSV RQQHKNFGLWGSLKNITWGN
Source	in vitro E.coli expression system
Target Names	tim21
Protein Names	Recommended name: Mitochondrial import inner membrane translocase subunit tim21
Expression Region	24-223
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.