



Recombinant Arabidopsis thaliana Reticulon-like protein B11 (RTNLB11)

Product Code	CSB-CF878578DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9LT71
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MGDSVSSSSRVSVHQSLGAGSVADLLLWRNRTGAVILLISSTGFWFLFERAGY NLLSFVS NVLLLLVAIFFLWAKSATVLNRPLPPVPNMEIPEEFANKAADDLRVWINYVLSIA SDITI ARNPIRLLQVSLVLWAI SYVGT LINS LTLVYIGVLLSLSFP IVYEKYQDHIDEKVNS TSK FVRNISRKFPLPINKEKKHQ
Source	in vitro E.coli expression system
Target Names	RTNLB11
Protein Names	Recommended name: Reticulon-like protein B11 Short name= AtRTNLB11
Expression Region	1-200
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.