



Recombinant *Oryza sativa* subsp. *indica* UPF0603 protein Osl_019212, chloroplastic (Osl_19898)

Product Code	CSB-CF386094OFF
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A2Y4G9
Product Type	Transmembrane Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>indica</i> (Rice)
Sequence	SEFDVLNNGGPPEDTYVVDDAGVLSRVTKSDVKRLVRDLESRKNIRINFITVRKL TSKADA FEYADQVLEKWYPTVEEGNNGIVVLVTSQKEGAI TGGPAFVQAVGDEILDST VSENLPV LATDEKYNEAIYTTAKRLAAAIDGLPDPGGPTFKDNKRESNFKTKEETEEKRGQ FTLVVG GLLVIAFVVPMAQYYAYISKK
Source	in vitro E.coli expression system
Target Names	Osl_19898
Protein Names	Recommended name: UPF0603 protein Osl_019212, chloroplastic
Expression Region	99-299
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.