



# Recombinant *Oryza sativa* subsp. *japonica* UPF0603 protein Os05g0401100, chloroplastic (Os05g0401100, LOC\_Os05g33280)

<b>Product Code</b>	CSB-CF738112OFG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q6ATY4
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Sequence</b>	SEFDVLNGGPPEDTYVDDAGVLSRVTKSDVKRLVRDLESRKNIRINFITVRKL TSKADA FEYADQVLEKWYPTVEEGNNKGIVVLVTSQKEGAI TGGPAFVQAVGDEILDST VSENLPV LATDEKYNEAIYTTAKRLAAAIDGLPDPGGPTFKDNKRESNFKTKEETEEKRGQ FTLVVG GLLVIAFVVPMAQYYAYISKK
<b>Source</b>	in vitro <i>E.coli</i> expression system
<b>Target Names</b>	Os05g0401100
<b>Protein Names</b>	Recommended name: UPF0603 protein Os05g0401100, chloroplastic
<b>Expression Region</b>	99-299
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.