



Recombinant Mouse Receptor expression-enhancing protein 1 (Reep1)

Product Code	CSB-CF803856MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BGH4
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MVSWIISRLVVLIFGTLYPAYYSYKAVKSKDIKEYVKWMMYWIIFALFTTAETFT DIFLC WFPFYELKIAFVAWLLSPYTKGSSLLYRK FVHPTLSSKEKEIDDCLVQAKDRS YDALVH FGKRGLNVAATAAVMAASKGQGALSERLRSFSMQDLTTIRGDGAPAPSGPPP PGTGRSSG KHSQPKMSRSASESAGSSGTA
Source	in vitro E.coli expression system
Target Names	Reep1
Protein Names	Recommended name: Receptor expression-enhancing protein 1
Expression Region	1-201
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.