



Recombinant Chicken E3 ubiquitin-protein ligase RNF152 (RNF152)

Product Code	CSB-CF019848CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	E1C2W7
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	METLSQDSLLECQICFNYYSPRRRPKLLDCKHTCCSVCLQQMRTSQKDLRCP WCRGITKL PPGYSVSQLPDDPEVIAVIAIPHTSEHTPVFIKLP SNGCYMLPLPLSKERAMLPG DIGCR LLPGSQQKSLAVVTIPAEQQPLQGGLPAEAGAEEDRRGVVKSSTWSGVCTVI LVACVLV FLLGIVLHNMSCISKRFTVISCG
Source	in vitro E.coli expression system
Target Names	RNF152
Protein Names	Recommended name: E3 ubiquitin-protein ligase RNF152 EC= 6.3.2.- Alternative name(s): RING finger protein 152
Expression Region	1-203
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.