



Recombinant Human Low-density lipoprotein receptor class A domain-containing protein 1 (LDLRAD1)

Product Code	CSB-CF729206HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5T700
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MNKVFPQGGENGYTAAESKAHPGGEAGGGHLCCSRRGACLSASLLLLLATVAA LIALVTIL GLPSCTPGAQACITLTNRTGFLCHDQRSCIPASGVCDGVRTCTHGEDEDESLC RDVPQSL PHFLVAHCGDPASWIYSDQKCDGTNNGDCSDELSPVTVCPGPGWWRRC PSTFFKYCDC IPRHLCDRHVQHCSDWSDEYACPGP
Source	in vitro E.coli expression system
Target Names	LDLRAD1
Protein Names	Recommended name: Low-density lipoprotein receptor class A domain-containing protein 1
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.