



Recombinant *Ashbya gossypii* Golgi to ER traffic protein 1 (GET1)

Product Code	CSB-CF762201DOT
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q754R6
Product Type	Transmembrane Protein
Immunogen Species	<i>Ashbya gossypii</i> (strain ATCC 10895 / CBS 109.51 / FGSC 9923 / NRRL Y-1056) (Yeast) (<i>Eremothecium gossypii</i>)
Sequence	MDYWILLVLAFLVADKSWHLTGLLATKLTSPERLQQILIRERQELHQQQQSLSA QDHYAKW TKNNRRDLVLDRIARVRKNYLESVEATKARLAKLKLKVVTVPFTALKFYKPKL PVYALP KGMFPRFIEGTLEHGWLYMALAPLNMKQFSEGASVAVSLGIWLFALLRVLGAI EFVLETLREQNPQVATETAKVHARTAQAASAN
Source	in vitro <i>E.coli</i> expression system
Target Names	GET1
Protein Names	Recommended name: Golgi to ER traffic protein 1 Alternative name(s): Guided entry of tail-anchored proteins 1
Expression Region	1-205
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.