



Recombinant Arabidopsis thaliana Cytochrome c oxidase assembly protein COX11, mitochondrial (COX11)

Product Code	CSB-CF814090DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8GWR0
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	STHSPSETKSQKMLYYLTAVVFGMVGLTYAAVPLYRTFCQATGYGGTVQRKE TVEEKIAR HSESGTVTEREIVVQFNADVADGMQWKFTPTQREVRVKPGESALAFYTAENK SSAPITGV STYNVTPMKAGVYFNKIQCFCFEEQRLLPGEQIDMPVFFYIDPEFETDPRMDGI NNLILS YTFFKVSEENTTETVNNNNSVPVQETN
Source	in vitro E.coli expression system
Target Names	COX11
Protein Names	Recommended name: Cytochrome c oxidase assembly protein COX11, mitochondrial
Expression Region	81-287
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.