



Recombinant Neisseria gonorrhoeae UPF0070 protein NGO0425 (NGO0425)

Product Code	CSB-CF528335NGA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O87406
Product Type	Transmembrane Protein
Immunogen Species	Neisseria gonorrhoeae (strain ATCC 700825 / FA 1090)
Sequence	MAAHLEEQQELDNFKYFWKTTGKWLFALLILAALGYLGTVYQNRAASQSQE AAAVLANI VEKAQNKAPQSEINAELSKLQQSYPHSISAAQATLMAAATEFDAQRYDVAEGH LKWVLSN QKDSLIALAAQRLGVVLLQQKKYDAALAALDTPVEADFAPLLMETKGDVYAA QEKSQEA LKNYGGQALEKMPQDSVGRELLQMKLDSLK
Source	in vitro E.coli expression system
Target Names	NGO0425
Protein Names	Recommended name: UPF0070 protein NGO0425 Alternative name(s): ORF2
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.