



Recombinant Bovine Probable glutathione peroxidase 8 (GPX8)

Product Code	CSB-CF649918BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2NL01
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MEPLTAYPLRCSGPKAKAFVLLSMVLCTVMLFLLQLKFLKPKINSFYTFEVKD ANGRVV SLEKFKGKVALVVNVASDCQLTDRNYLALQELHKEFGPFHFSVLAFPCNQFGE SEPRPSK EVVSFARNNFGVTFPIFHKIKILGSEAEPAFRFLVDSSKKEPRWNFWKYLVNPE GQVVKT WRPEEPIEVIRPEIAALIRQMIKKKEDL
Source	in vitro E.coli expression system
Target Names	GPX8
Protein Names	Recommended name: Probable glutathione peroxidase 8 Short name= GPx-8 Short name= GSHPx-8 EC= 1.11.1.9
Expression Region	1-209
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.