



Recombinant Arabidopsis thaliana Putative RING-H2 finger protein ATL50 (ATL50)

Product Code	CSB-CF875451DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9FHG8
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MDQKSDSFLSVSSISFSYSSSTDKDFDLICMISPIVLLYITLLSIIFVVAALIHLLVK FL HRPQTRLDDAYDGITESSTALQGRYQTRFNLHDAEIDQSFIDALPLLHYKTMIG LRHDL DCAVCLREFTAEDLRLLPKCSHAFHVECIDTWLLTNSTCPLCRDNLLLLGLTG TASSST IVLVHESDGDNSQSDSSFMLTDLDDVESK
Source	in vitro E.coli expression system
Target Names	ATL50
Protein Names	Recommended name: Putative RING-H2 finger protein ATL50
Expression Region	1-210
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.