



Recombinant Vaccinia virus Protein F9 (F9)

Product Code	CSB-CF303890VAD
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P68453
Product Type	Transmembrane Protein
Immunogen Species	Vaccinia virus (strain L-IVP) (VACV)
Sequence	MAETKEFKTLYNLFIDSYLQKLAQHSIPTNVTCIAIHIGEIVIGQFKNCALRITNKCM SNSR LSFTLMVESFIEVISLLPEKDRRAIAEEIGIDLDDVPSAVSKLEKNCNAYAENVNII DIQ KLDIGECSAPPGQHMLLQIVNTGSAEANCGLQTIVKSLNKIYVPPIIENRLPYD PWFLV GVAILLVIFTVAICSI RRNLALKYRYGTFLYV
Source	in vitro E.coli expression system
Target Names	F9
Protein Names	Recommended name: Protein F9
Expression Region	1-212
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.