



Recombinant *Saccharomyces cerevisiae* Altered inheritance of mitochondria protein 36, mitochondrial (AIM36)

Product Code	CSB-CF513185SVN
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	C8ZF15
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain Lalvin EC1118 / Prise de mousse) (Baker's yeast)
Sequence	SSTDSSTKRSNKSDKIDAPGFKKIFLVAIIGTVIFVKTQSLDKNKPKTTLSEEEF ENVV KGLKRRVAIFPQGEVDIKFSLSPSIEETRKVLQKSQGDDINELQFVDPVKVIDYY RTLTD DRYEALLNEYKYYGCDTYAYNLPTGMLVMLLGRYFKENFKAGDKLVVVNFP HSIADATR FENEVSIVSKIFVPRKLSGSDVCKYYETVVKADII
Source	in vitro E.coli expression system
Target Names	AIM36
Protein Names	Recommended name: Altered inheritance of mitochondria protein 36, mitochondrial Alternative name(s): Found in mitochondria protein 39
Expression Region	41-255
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.