



Recombinant Arabidopsis thaliana Uncharacterized mitochondrial protein AtMg00910 (AtMg00910)

Product Code	CSB-CF309365DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92528
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MPTANQLIRHGREEKRRTDRTEVLVFGLLVTRIIRFVHSVLFPIPVFCSIKVLDDY FCSL PIIDKLSKKWQLIWFYVLSVILCKSLFAVGYLWMDDLRAISQFYPVVSGLGG GNTMP PTNPSEGGLEGGYAHENEHSHDQQRGSPFWSKEYKESGSKRFLNLEVED QNTDTIGE QVKAESGKCEKIKAKIIAKTHELLVSEDTKFQIKTI
Source	in vitro E.coli expression system
Target Names	AtMg00910
Protein Names	Recommended name: Uncharacterized mitochondrial protein AtMg00910 Alternative name(s): ORF215a
Expression Region	1-215
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.