



Recombinant Uncharacterized protein yidX (yidX)

Product Code	CSB-CF359973SZB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P0ADM7
Product Type	Transmembrane Protein
Immunogen Species	Shigella flexneri
Sequence	MKLNFKGFFKAAGLFPLALMLSGCISYALVSHTAKGSSGKYQSQSDTITGLSQ AKDSNGT KGYVFGESLDYLITDGADDIVKMLNDPALNRHNIQVADDARFVLNAGKKKFT GTISLYY YWNNEEEKALATHYGFACGVQHCTRSLENLKGTIHEKNKNMDYSKVMFYHP FKVRFY EY YSPRGIPDGVSAALLPVTVTLDIITAPLQFLVYAVNQ
Source	in vitro E.coli expression system
Target Names	yidX
Protein Names	Recommended name: Uncharacterized protein yidX
Expression Region	1-218
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.