



Recombinant Bovine BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like (BNIP3L)

Product Code	CSB-CF666194BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3T013
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MSSHLVEQPPPPHNNNNNCEEQEQLPPPAGLNSSWVELPMNSSNGNDNG NGKNGGLEHV PSSSSIHNGDMEKILLDAQHESGQSSSRGSSHCDSPPSPQEDGQIMFDVEMHT SKDHSSQS EEEVAEGEKEVDALKKSVDWVSDWSSRPENIPPKEFHFRHPKRSVLSMRKS GAMKKGGI FSAEFLKVFIPSLFLSHVLALGLGIYIGKRLSTPSASTY
Source	in vitro E.coli expression system
Target Names	BNIP3L
Protein Names	Recommended name: BCL2/adenovirus E1B 19 kDa protein-interacting protein 3-like
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.