



Recombinant Human V-set and transmembrane domain-containing protein 1 (VSTM1)

Product Code	CSB-CF747633HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6UX27
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	YEDEKKNEKPPKPSLHAWPSSVVEAESNVTLKCAHSQNVTFVLRKVNDSGY KQEQSSAE NEAEFPFTDLKPKDAGRYFCAYKTTASHEWSESSEHLQLVVTDKHDELEAPS MKTDRIT FVAIFSCISILLFLSVFIIYRCSQHSSSSSEESTKRTSHSKLPEQEAAEADLSNME RVSL STADPQGVTY AELSTSALSEAASDTTQEPGGSHEYAALKV
Source	in vitro E.coli expression system
Target Names	VSTM1
Protein Names	Recommended name: V-set and transmembrane domain-containing protein 1
Expression Region	17-236
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.