



Recombinant *Drosophila melanogaster* Probable mitochondrial import inner membrane translocase subunit Tim17 3 (Tim17a1)

Product Code	CSB-CF895136DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9VGA2
Product Type	Transmembrane Protein
Immunogen Species	<i>Drosophila melanogaster</i> (Fruit fly)
Sequence	MEYNRQPCPIRIVEDCGCAFMMGTIGGSLFEFLKGFRNAPTGLQRRRLYGGIDL VKMRTPS IAGSFAVWGATFSTVDCALVHYRQREDAWNSILSGAATGGILAARNGIRAMAN SALVGCL VLAMIEGAGAAVATINAADKGAGIVIKPQRAQWEAILETIDPKRASSTQDFALAE FERVL DKCRASREPNULLQDIPVKSHERDSKQKPFYSLLDLVKLSQMF
Source	in vitro <i>E.coli</i> expression system
Target Names	Tim17a1
Protein Names	Recommended name: Probable mitochondrial import inner membrane translocase subunit Tim17 3
Expression Region	1-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.