



Recombinant Uncharacterized protein ZK512.1 (ZK512.1)

Product Code	CSB-CF333919CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34639
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	ATWVKITKPTITWEDMEVSLEGEKKFFCGEFDSSFTGRYHWRFNGSSILPERT QIHRNQF VFLAGANAIRNQLPGEYECCVRETLGNACYSRMIVVQNRTDNHNIDMTNSTLL LADEGNT YYIRMHDVKRVEGVKCTLDGVNVDNFKYPFLGRQTKKTPYHLKIENIERGGE VNCDLRL HKKEIVQKTFDIRLLRGFISSQLPQFVYLIVFTIIGYILRL
Source	in vitro E.coli expression system
Target Names	ZK512.1
Protein Names	Recommended name: Uncharacterized protein ZK512.1
Expression Region	19-240
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.