



Recombinant Danio rerio TM2 domain-containing protein 3 (tm2d3)

Product Code	CSB-CF023614DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A5PLF5
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	YLSSPHVGQDPPYLAQQKSVMSPPVALTASSASPVVTDNYVSKCPSGGLCSK LPADCIIC ALHHNCSYGRPHNYTCRPRAGVHCVSDQGERQQNFTLSLLCRFCFQLDASQ YRCSNSSDC MTVSCPRRRRYNASCEVLEHVHCLGKRVFQKRLFCNWTGGYKWSTALALSITL GGFGADRF YLGQWREGLGKLFSGGLGIWTLIDVLLIGVGYVGPADGSLYI
Source	in vitro E.coli expression system
Target Names	tm2d3
Protein Names	Recommended name: TM2 domain-containing protein 3
Expression Region	22-244
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.