



Recombinant *Saccharomyces cerevisiae* Altered inheritance rate of mitochondria protein 38 (AIM38)

Product Code	CSB-CF345699SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P53721
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Sequence	MKILTQDEIEAHRSHLTKGGIEGALAGFAISAIIFKVLPRRYPKFKPSTLTWSIKT ALWI TPPTVLTAICAEAEASNNFDATMYGSGSSSEDALDEHRRWKSLSLTKDKFVEGLS NNKYKII TGAWAASLYGSWVIVNKDPIMTKAQKIVQARMYAQFITVGLLLASVGLSMYEN KLHPNKQ KVNEMRRWENALRVAEEEEERLEKEGRRTGYVSNEERINSKIFKS
Source	in vitro E.coli expression system
Target Names	RCF2
Protein Names	Recommended name: Altered inheritance rate of mitochondria protein 38
Expression Region	1-224
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.