



Recombinant Human coronavirus NL63 Non-structural protein 3 (3)

Product Code	CSB-CF754294HIX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6Q1S1
Product Type	Transmembrane Protein
Immunogen Species	Human coronavirus NL63 (HCoV-NL63)
Sequence	MPFGGLFQLTLESTINKSVANLKLPPHDVTVLRDNLKPVTTTSTITAYLLVSLFV TYFAL FKPLTARGRVACFVLKLLTLFVYVPLLVLFGMYLDSFIIFSTLLFRFIHVGYYAYL YKNF SFVLFNVTKLCFVSGKWCWYLEQSFYENRFAAIYGGDHYVVLGGETITFVSFDDL YVAIRG SCEKNLQLMRKVDLYNGAVIYIFAEEPVVGIVYSSQLYEDVPSIN
Source	in vitro E.coli expression system
Target Names	3
Protein Names	Recommended name: Non-structural protein 3 Short name= Protein 3 Short name= ns3 Alternative name(s): Accessory protein 3a
Expression Region	1-225
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.