



Recombinant Danio rerio Leucine-rich repeat-containing protein 3B (Irrc3b)

Product Code	CSB-CF013137DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	A3KNN3
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	CPKGCTCQRSESPPHGLNVTCSLSRLKEIPPDVPPDTQLLQLDRNHISLVPDRI FHGLRM LRRNLNSHNAVETLGEGAFIGLEGSLEVLDLSYNRITSVHKDAFARLKARVVVD NNPWHC DCALQQALGGMAHNHERVLCRSSELRDQEGQPFAVDADLCNLAKRTTDYA MLVTMFGWF AMVISYVVVYVRQNQEDARRHLEYLKSLPSKPKKPDEPEDISTVV
Source	in vitro E.coli expression system
Target Names	Irrc3b
Protein Names	Recommended name: Leucine-rich repeat-containing protein 3B
Expression Region	34-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.