



Recombinant Bovine ADP-ribosylation factor-like protein 6-interacting protein 6 (ARL6IP6)

| | |
|--------------------------|--|
| Product Code | CSB-CF659524BO |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | Q3MHM8 |
| Product Type | Transmembrane Protein |
| Immunogen Species | Bos taurus (Bovine) |
| Sequence | MSFVESGRRSAPLRRRPGTPVPFARPAYSVFSQGDSWGEGEVEEEEVGGCDQV ARDLRAEFS AGSSSKLRKDPVLQPDGDGSPVLPDKRNGIFSADAGGKALARRWPVQVLSIL CSLLFAIL LACLLAITYLIVKELHAENLKNEDDVNTGLLGFWSLLIISLTAGFSCCSFSWTVTY FDSF EPGMFPPTPLSPARFKKMTGHSFHMGYSMAILNGIVAALTVAWCLM |
| Source | in vitro E.coli expression system |
| Target Names | ARL6IP6 |
| Protein Names | Recommended name: ADP-ribosylation factor-like protein 6-interacting protein 6 Short name= ARL-6-interacting protein 6 Short name= Aip-6 |
| Expression Region | 1-226 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |