



Recombinant Mouse ADP-ribosylation factor-like protein 6-interacting protein 6 (Arl6ip6)

Product Code	CSB-CF805886MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BH07
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MSFVESWRFAGARRRRQVTPGPATRPGYSDYTQGDSWGE GEGDEDEGCD QVARDLRAEFS ARASSETKRAPLLPRVGDGSPVLPDKRNGIFPATAAKRTQARRWPIQALSILCS LLFAVL LAFLLAIAYMIVKELHAENLKNEDDIHTGLLGFWSLLIISLTAGLSCCSFSWTVTY FDSF EPGMFPPTPLSPARFKKLTGHSFHMGYSMAILNGIVAALTVAWCLM
Source	in vitro E.coli expression system
Target Names	Arl6ip6
Protein Names	Recommended name: ADP-ribosylation factor-like protein 6-interacting protein 6 Short name= ARL-6-interacting protein 6 Short name= Aip-6
Expression Region	1-226
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.