



Recombinant Human Leucine-rich repeat-containing protein 3B (LRRC3B)

Product Code	CSB-CF857000HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q96PB8
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	CPKGCLCSSSGGLNVTCSNANLKEIPRDLPPETVLLYLDSNQITSIPNEIFKDLH QLRVL NLSKNGIEFIDEHAFKGV AETLQTLDSLSDNR IQSVHKNAFN NLKARARIANNPW HCDCTL QQVLRSMASNHETAHN VICKTSVLDEHAGRPFLNAANDADLCNL PKKTTDYA MLVTMFGW FTMVISYVVYYVRQNQEDARRHLEYL KSLPSRQKKADEPDDISTVV
Source	in vitro E.coli expression system
Target Names	LRRC3B
Protein Names	Recommended name: Leucine-rich repeat-containing protein 3B Alternative name(s): Leucine-rich repeat protein LRP15
Expression Region	34-259
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.