



# Recombinant *Malacothrix typica* Cytochrome c oxidase subunit 2 (MT-CO2)

<b>Product Code</b>	CSB-CF657792MCAN
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q37547
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Malacothrix typica</i> (Long-eared mouse) (Gerbil mouse)
<b>Sequence</b>	MAYPLQLGLQDATSPIMEELMNFHDHTLMIVFLISSLVLYVISSMLTTKLHTST MDAQE VETIWTILPAVILIMIALPSLRILYMMDEINNPVLTVKTMGHQWYWSYEYTDYED LCFDS YMIPTTDLKPGEFRLLEVDNRVILPMELPIRMLISSEDLHLSWAIPSLGLKTDIAIP GRLN QATISSNRPGLFYGQCSEICGSNHSFMPPIILEMVPLKNFETWSVSMI
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	MT-CO2
<b>Protein Names</b>	Recommended name: Cytochrome c oxidase subunit 2 Alternative name(s): Cytochrome c oxidase polypeptide II
<b>Expression Region</b>	1-227
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.