



Recombinant Arabidopsis thaliana Peroxisomal membrane protein 11B (PEX11B)

Product Code	CSB-CF871163DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9STY0
Product Type	Transmembrane Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Sequence	MSLDTVDKLVVFLAKRDGIDKLVKTFQYVAKLACWHVEATRPEAADRFKKWEV ASGLSRK AFRTGRSLTGFNALRRNPGATPMIRFLAVLANSGEMVYFFFDHFLWLSRIGSID AKLAKK MSFISAFGESFGYTFFIIDCIFIKQRLKSLKKLQHSTDEPKKEEIGAKISEIRGDIVM RL MGISANVADLLIALAEIHPNPFNHTITLGISGLVSAWAGWYRNWPS
Source	in vitro E.coli expression system
Target Names	PEX11B
Protein Names	Recommended name: Peroxisomal membrane protein 11B Alternative name(s): Peroxin-11B Short name= AtPEX11b
Expression Region	1-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.