



Recombinant Human Probable N-acetyltransferase 8B (NAT8B)

Product Code	CSB-CF883394HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9UHF3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MAPYHIRKYQESDRKSVVGLLSGGMAEHAPATFRLLKLPRTLILLGGALALL LVSGSW ILALVFSLSLLPALWFLAKKPWTRYVDIALRTDMSDITKSYLSECGSCFWVAES EEKVVG TVGALPVDDPTLREKRLQLFHLSVDNEHRGQGIKALVRTVLQFARDQGYSEV VLDTSTNI QLSAMGLYQSLGFKKTGQSFFHVWARLVDLHTVHFIYHLPQAQAGRL
Source	in vitro E.coli expression system
Target Names	NAT8B
Protein Names	Recommended name: Probable N-acetyltransferase 8B EC= 2.3.1.- Alternative name(s): Camello-like protein 2
Expression Region	1-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.