



Recombinant *Macaca fascicularis* Cytochrome c oxidase subunit 2 (MT-CO2)

Product Code	CSB-CF015073MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P11948
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	MAHPVQLSLQDATSPIMEELITFHDHAFMAMSLISFLVLYALALTLTTKLTNTNIT DAQE METIWTILPAVILILIALPSLRVLYLTDEVNDPSLTIKSIGHQWYWTYEYTDYGGLI FNS YMLPPLFLKPGDLRLLEVDNRVVLPIEAPVRMMITSQDVLHSWTIPTLGLKTD VPGRLN QTVFTATRPGVYYGQCSEICGANHSFMPIVADLIPLKIFEMGPVFTL
Source	in vitro <i>E.coli</i> expression system
Target Names	MT-CO2
Protein Names	Recommended name: Cytochrome c oxidase subunit 2 Alternative name(s): Cytochrome c oxidase polypeptide II
Expression Region	1-227
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.