



Recombinant Danio rerio Ankyrin repeat domain-containing protein 46 (ankrd46)

Product Code	CSB-CF757762DIL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6TNT2
Product Type	Transmembrane Protein
Immunogen Species	Danio rerio (Zebrafish) (Brachydanio rerio)
Sequence	MSYVFINDSSQTSVPLLQACIDGDL SFARRLLETGCDPNIRDHRGRTGLHLAAA RGNVDI CRFLHKFGADLLATDYQGNTALHLCGHVDTIQFLVSNGLKIDICNHNGSTPLVL AKRRGV NKDAIRLLEGLEEQEVKGFNRGAHASKLEAMQMAESESAMESHLLNPNLQNS EGVLSSFR STWQEFVEDLGFWRVLLLLLVIALLSLGIAYYVSGVLPFSASQLELVH
Source	in vitro E.coli expression system
Target Names	ankrd46
Protein Names	Recommended name: Ankyrin repeat domain-containing protein 46
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.