



Recombinant *Adalia bipunctata* Cytochrome c oxidase subunit 2 (COII)

Product Code	CSB-CF015073ACZ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29871
Product Type	Transmembrane Protein
Immunogen Species	<i>Adalia bipunctata</i> (Two-spotted ladybird beetle) (<i>Coccinella bipunctata</i>)
Sequence	MSTWKSSLFLDSSSFLLQLRFFHHDHALLILNMITGAVAYIMISLLFNKYNHRFLL EGHT VETIWTILPAFTLIFIALPSLKLILYIDEIRNPLVTLKTIGHQWYWTYEYSDFKKLEF DS YMLSVDNLNPFNFRLLLEVDNRTILPYLSNIRLLTSSADVIHSWTIPSSGVKIDASP GRLN QMSFTLNRTGIFYGQCSEICGANHSFMPISIESISPSNFIKWVNKSSL
Source	in vitro E.coli expression system
Target Names	COII
Protein Names	Recommended name: Cytochrome c oxidase subunit 2 EC= 1.9.3.1 Alternative name(s): Cytochrome c oxidase polypeptide II
Expression Region	1-228
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.