



Recombinant Escherichia coli Protein gltF (gltF)

Product Code	CSB-CF333377ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P28721
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	TDSTDTELTIIIGEYTPGACTPVVTGGGIVDYGKHHNSALNPTGKSNKLVQLGRK NSTLNI TCTAPTLIAVTSKDNRQSTIVALNDTSYIEKAYDTLVDMKGTKNAFGLGSAPNG QKIGAA SIGIDRSNGGIHAADDTGEIPVDLIQTDHWSAATPTWKASSNGAFCSLTSCSAI ERGYSV AKTGELTPVAITAVTFPLLIDAAVNDNTILGSDETILKLDGNVTISVQYL
Source	in vitro E.coli expression system
Target Names	gltF
Protein Names	Recommended name: Protein gltF
Expression Region	26-254
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.