



Recombinant Toxoplasma gondii Dense granule protein 6 (GRA6)

Product Code	CSB-CF631656TOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q27003
Product Type	Transmembrane Protein
Immunogen Species	Toxoplasma gondii
Sequence	MAHGGIHLRQKRNFVTVSTVAVVFVFMGVLVNSLGGVAVAADSGGVKQT PSETGSSG GQQEAVGTTEDYVNSSAMGGGQGDLSLAEDDTTSEAAEGDVDPFPVLANEGK SEARGPSLE ERIEEQGTRRRYSSVQEPQAKVPSKRTQKRHRLIGAVVLAVSVAMLTAFFLRR TGRRSPQ EPSGDGGGNDAGNNAGNGGNEGRGYGGRGEGGAEDDRRPLHPERVNVFD Y
Source	in vitro E.coli expression system
Target Names	GRA6
Protein Names	Recommended name: Dense granule protein 6 Short name= Protein GRA 6Alternative name(s): Antigen p32 Protein p33
Expression Region	1-230
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.