



Recombinant Escherichia coli UPF0702 transmembrane protein ycaP (ycaP)

Product Code	CSB-CF303456ENV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P75839
Product Type	Transmembrane Protein
Immunogen Species	Escherichia coli (strain K12)
Sequence	MKAFDLHRMAFDKVPFDLGEVALRSLYTFVLVFLFLKMTGRRGVRQMSLFEV LIILTLG SAAGDVAFYDDVPMVPVLIVFITLALLYRLVMWLMHSEKLEDLLEGKPVVIED GELAW SKLNNSNMTEFEFFMELRLRGVEQLGQVRLAILETNGQISVYFFEDDKVKPGLL ILPSDC TQRYKVPESADYACIRCSEIIHMKAGEKQLCPRCANPEWTKASRAKRV
Source	in vitro E.coli expression system
Target Names	ycaP
Protein Names	Recommended name: UPF0702 transmembrane protein ycaP
Expression Region	1-230
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.