



Recombinant Zea mays Chlorophyll a-b binding protein, chloroplastic

Product Code	CSB-CF356910ZAX
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P06671
Product Type	Transmembrane Protein
Immunogen Species	Zea mays (Maize)
Sequence	RKTAAKAKPAAASGSPWYGPDRVLYLGPLSGEPPSYLTGEFPGDYGWDTAG LSADPETFA KNRELEVIHCRWAMLGALGCVFPELLARNGVKFGEAVWFKAGSQIFSEGGLD YLGNP SLI HAQSILAIWACQVVL MGAVEGYRIAGG PLGEVVDPLYPGGSFDPLGLADDPEA FGELKVK ELKKGRLAMLSMFGFFVQAIVTGKGP L ENLADHIADPVNNNAWAYATNFVPGK
Source	in vitro E.coli expression system
Protein Names	Recommended name: Chlorophyll a-b binding protein, chloroplastic Alternative name(s): LHCII type I CAB Short name= LHCP
Expression Region	33-265
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.